


# Innovations in Genomic Analysis:

## Downstream analysis of Illumina Sequencing Data

Marco Cappelletti  
Product Marketing Manager  
Europe

# Agenda


- Genome Analyzer IIx Sequencing Technology
- Applications overview
- GA/IIx SW Improvements
- Downstream Data Analysis




# The Genome Analyzer<sub>IIx</sub> and Software Advancements

*65% Increase in data output*

- 20-25 Gb of high quality data / run
- 2.5 Gb / day
- >300M reads per paired-end run
- 2 x 75bp supported read length
- Raw Accuracy:  $\geq 98.5\%$
- Short-Insert PE and long-insert Mate Pair sequencing



# Illumina Genome Analyzer: A paradigm shift



# Simplest Sequencing Process - 5 to 10 days WF

## 1 *Library prep (~ 6 hrs)*



Fragment DNA  
↓  
Repair ends / Add A overhang  
↓  
Ligate adapters  
↓  
Select ligated DNA

## 2 *Automated Cluster Generation (~ 5 hrs)*



Up to 96 samples

Hybridize to flow cell  
↓  
Extend hybridized oligos  
↓  
Perform bridge amplification


## 3 *Sequencing (~ 4-9 days\*)*



Up to 96 samples

Perform sequencing on forward strand  
↓  
Re-generate reverse strand  
↓  
Perform sequencing on reverse strand

# Flow cell







# Genomic DNA Library Prep

1

DNA fragments



PCR  
6 – 15 cycles





Library

# Cluster Generation

## *Cluster station*

- Aspirates DNA samples into flow cell
- Automated amplified clonal clusters






# Cluster Generation

## *Hybridize Fragment & Extend*


- >150 M single molecules hybridize to the lawn of primers
- Bound molecules are then extended by polymerases



# Cluster Generation


## *Denature Double-stranded DNA*

- Double-stranded molecule is denatured
- Original template is washed away
- Newly synthesized covalently attached to the flow cell surface



# Cluster Generation


## *Covalently-Bound Spatially Separated Single Molecules*



# Cluster Generation

## *Bridge Amplification*


- Single-strand flips over to hybridize to adjacent primers to form a bridge
- Hybridized primer is extended by polymerases



# Cluster Generation

## *Bridge Amplification*


- Double-stranded bridge is formed



# Cluster Generation

## *Bridge Amplification*

- Double-stranded bridge is denatured
- Result: Two copies of covalently bound single-stranded templates






# Cluster Generation

## *Bridge Amplification*


- Single-strands flip over to hybridize to adjacent primers to form bridges
- Hybridized primer is extended by polymerase



# Cluster Generation


## *Bridge Amplification*

- Bridge amplification cycle repeated until multiple bridges are formed



# Cluster Generation


- dsDNA bridges denatured
- Reverse strands cleaved and washed away






# Cluster Generation

- Leaving a cluster with forward strands only




# Genome Analyzer Sequencing reaction

- Sequencing primer is hybridized to adapter sequence




# Sequencing






# Sequencing



~1000 copies  
per cluster



**>150 Million Clusters Per Flow Cell**

# Broadest range of applications

*Optimized, streamlined and easy-to-use reagent solutions*

## Sample Prep

### Whole genome

- Resequencing
- De-novo
- Targeted
- Metagenomics


### Transcriptome

- RNA-Seq
- DGE
- Small RNA
- miRNA

### Regulation

- Methylation
- ChIP-Seq

## Automated Cluster Generation



## Sequencing

**Epigenomics**

**Transcriptomics**

***de novo* sequencing**

**Metagenomics**

**Genomics**

**...go where the biology takes you**

# Structural Variation Analysis

Epigenomics


Transcriptomics


*de novo* sequencing

Metagenomics


Genomics

...to go where the biology takes you





# Reference



Deletion Inversion! Deletion

Local *de novo* assembly

Long Inserts

Deletion?

Short Inserts

Inversions?

Middle

Normal?

**Epigenomics**

**Transcriptomics**

***de novo* sequencing**

**Metagenomics**

**Genomics**



# Nature Publications

Author	Number of Reads	Average Length
Shi, <i>et al.</i>	389K	97 bp
Dinsdale, <i>et al.</i>	14.6M	105 bp
Mou, <i>et al.</i>	307K	96 bp
Warnecke, <i>et al.</i>	300K	100 bp
Turnbaugh, <i>et al.</i>	1.7M	93 bp

# Nature Publications

Paper	Number of Reads	Average Length
Metatranscriptomics	389K	97 bp
Biome Profiling	14.6M	105 bp
Ocean Bacteria	307K	96 bp
Termite Microbiota	300K	100 bp
Obesity Microbiome	1.7M	93 bp

# Nature Publications

**Reads per run:  
150-300 Million**

**Number  
of Reads**

**% of a 150M  
read run**



389K

0.3%

14.6M

9.7%

307K

0.2%

300K


0.2%

1.7M

1.1%

# The Value of Long, Overlapping Reads

Error Rate



**Epigenomics**

**Transcriptomics**

***de novo* sequencing**

**Metagenomics**

**Genomics**


“

When the read length exceeds a certain threshold, the ***read length barrier***, the efficiency reaches nearly 100%, so that the read length indeed does not matter.

”

Chaisson et al., De novo *Fragment Assembly with Short Mate-Paired Reads: Does the Read Length Matter?*, *Genome Research*, October 22, 2008







October, 2008



*Salmonella seftenberg*

Illumina

N50 contig size:	139,353
Largest contig:	395,600
Average contig:	63,969
Total bases of contigs:	4.80Mb
Coverage of genome:	99.8%



**Beijing Genome Institute**

January 9, 2009

- 3GB genome
- Paired 75-base reads
- >95% gene regions
- N50 contig: 300Kb

**Epigenomics**

**Transcriptomics**


***de novo* sequencing**


**Metagenomics**

**Genomics**




- Capabilities
- Length of read
- Short-insert Paired Ends
- Raw Read Accuracy






Source: Illumina R&D

## Ensembl



## Illumina GA






**Epigenomics**

**Transcriptomics**

*de novo* sequencing



**Metagenomics**

**Genomics**



# GA Applications Published 24 months from launch

<b>nature</b>	<b>42</b>
<b>Science</b>	<b>14</b>





# Genome Analyzer<sub>IIx</sub>

## Software Advancements


*Increased Output, Simplified Computing*

# Increased output with reduced computing infrastructure

*More gigabases of data for fewer gigabytes of computing power!*

## Sequencing Control Software v2.4


- Includes new Real Time Analysis (RTA) feature
- Image extraction and real time base calling on instrument computer
- Shorter time to results
  - Performed simultaneously with sequencing
  - Eliminates need to transfer images and intensities across network
  - Base calls and quality scores within hours of end of run



# New Software Delivers Up to 40% More Data Per Run


## ● Pipeline 1.4 – Enhanced analysis algorithm

- Increases yield, improves accuracy
- Improved cluster delineation
- More clusters pass filter
- Lower error rates



## ● SCS 2.4 – Integrated autofocus

- Easier to use, removes user error
- Less variability in focus quality




# What is *Real Time Analysis*?

- The RTA module analyzes data as it leaves the Genome Analyzer
  - **Produces base calls, including Phred-like quality scores**
  - **Generates reports, to assess run and library quality**
  - Performs image analysis, generation of cluster intensities
- RTA simplifies the data management process
  - **Eliminates the need to transfer images from computer to computer**
  - Includes optional mechanisms for complete or selective archiving of images
  - Includes optional mechanisms for archiving of intensities
- RTA improves the system performance
  - **Minimizes time to results – base calls and qualities generated within hours of the end of the run**
  - Removes dependencies on network availability
  - Minimizes the time spent analyzing data after the run




# Reversing the Trend

*Simplified computing, Smaller storage needs, Faster analysis*



# FireCrest module: Start From the Spots...




tiff image files

- Clusters identification and assigns intensities to them
- The output is a simple text file

s\_1\_0007\_int.txt - WordPad

File	Edit	View	Insert	Format	Help						
1	7	135	563	168.9	347.7	739.1	24966.8	202.2	299.7	207.0	21984.4
1	7	180	621	231.5	341.9	457.7	21423.8	229.3	382.9	16319.2	20217.5
1	7	245	626	218.4	356.8	501.6	21362.3	165.5	319.7	467.9	19749.5
1	7	241	509	187.7	382.7	537.4	20767.7	1489.2	10304.1	161.0	482.7
1	7	214	595	173.5	372.1	686.1	20302.4	8387.1	12746.0	158.4	540.8
1	7	155	544	172.2	339.5	538.3	19608.9	307.6	418.8	364.9	17172.9
1	7	301	507	353.8	672.1	782.0	26448.1	1881.2	12332.1	191.9	743.0
1	7	175	606	210.4	333.4	523.2	19248.3	164.4	308.7	535.9	20587.5
1	7	242	522	267.9	513.0	606.8	19056.7	6265.6	10442.1	1884.7	2440.8
1	7	196	522	220.2	455.9	486.6	18895.4	189.5	352.8	12299.4	14331.7
1	7	237	612	167.0	457.7	531.0	18835.2	713.8	992.0	416.4	18774.3
1	7	160	528	172.6	400.7	651.9	18686.9	1265.7	8500.6	241.3	524.1
1	7	164	543	205.7	385.0	488.4	18480.5	1410.3	9968.3	76.7	-343.0
1	7	179	581	207.2	372.9	560.1	18462.2	140.7	282.9	314.4	16462.8
1	7	226	623	218.3	400.6	474.6	18392.9	7333.1	10759.6	158.2	640.2
1	7	139	583	241.0	358.9	563.7	18183.9	226.9	302.0	11925.1	15357.5
1	7	220	618	223.1	496.8	553.2	18176.5	1338.5	10208.8	315.3	594.6
1	7	360	507	194.0	339.0	660.3	24628.4	294.7	590.6	620.8	26846.9
1	7	334	512	249.8	590.6	638.9	24101.4	6787.9	11276.9	602.5	177.3
1	7	155	517	218.7	345.4	554.6	17715.4	1415.3	8446.5	177.4	523.2
1	7	343	541	183.5	375.9	678.6	23803.5	6715.9	11488.7	189.9	684.9
1	7	241	608	208.6	361.2	457.0	17245.5	6250.2	9519.9	112.1	34.4
1	7	176	520	226.3	338.6	457.9	17172.1	179.5	300.5	387.3	16274.9
1	7	371	592	298.6	566.4	626.1	23249.9	6698.6	10982.2	146.3	210.1
1	7	271	508	175.8	391.5	567.5	23181.2	1502.2	11095.5	158.9	605.8
1	7	195	503	236.4	389.5	485.4	16827.3	6096.1	8300.3	189.5	5778.0
1	7	301	592	181.8	378.8	553.6	22568.7	8013.1	13222.2	889.6	1211.8
1	7	248	548	197.7	525.1	543.6	16512.2	1560.8	10651.3	175.3	508.9
1	7	242	532	208.7	386.0	508.1	16468.5	1355.0	8800.2	155.7	283.8



intensity files

# From Intensities to Reads

## Intensity Files

s\_1\_0007\_int.txt - WordPad

1	7	135	563	168.9	347.7	739.1	24966.8	202.2	299.7	207.0	21984.4
1	7	180	621	231.5	341.9	457.7	21423.8	229.3	382.9	16319.2	20217.5
1	7	245	626	218.4	356.8	501.6	21362.3	165.5	319.7	467.9	19749.5
1	7	241	509	187.7	382.7	537.4	20767.7	1489.2	10304.1	161.0	482.7
1	7	214	595	173.5	372.1	686.1	20302.4	8387.1	12746.0	158.4	540.8
1	7	155	544	172.2	339.5	538.3	19608.9	307.6	418.8	364.9	17172.9
1	7	301	507	353.8	672.1	782.0	26448.1	1881.2	12332.1	191.9	743.0
1	7	175	606	210.4	333.4	523.2	19248.3	164.4	308.7	535.9	20587.5
1	7	242	522	267.9	513.0	606.8	19056.7	6265.6	10442.1	1884.7	2440.8
1	7	196	522	220.2	455.9	486.6	18895.4	189.5	352.8	12299.4	14331.7
1	7	237	612	167.0	457.7	531.0	18835.2	713.8	992.0	416.4	18774.3
1	7	160	528	172.6	400.7	651.9	18686.9	1265.7	8500.6	241.3	524.1
1	7	164	543	205.7	385.0	488.4	18480.5	1410.3	9968.3	76.7	-343.0
1	7	179	581	207.2	372.9	560.1	18462.2	140.7	282.9	314.4	16462.8
1	7	226	623	218.3	400.6	474.6	18392.9	7333.1	10759.6	158.2	640.2
1	7	139	583	241.0	358.9	563.7	18183.9	226.9	302.0	11925.1	15357.5
1	7	220	618	223.1	496.8	553.2	18176.5	1338.5	10208.8	315.3	594.6
1	7	360	507	194.0	339.0	660.3	24628.4	294.7	590.6	620.8	26846.9
1	7	334	512	249.8	590.6	638.9	24101.4	6787.9	11276.9	602.5	177.3
1	7	155	517	218.7	345.4	554.6	17715.4	1415.3	8446.5	177.4	523.2
1	7	343	541	183.5	375.9	678.6	23803.5	6715.9	11488.7	189.9	684.9
1	7	241	608	208.6	361.2	457.0	17245.5	6250.2	9519.9	112.1	34.4
1	7	176	520	226.3	338.6	457.9	17172.1	179.5	300.5	387.3	16274.9
1	7	371	592	298.6	566.4	626.1	23249.9	6698.6	10982.2	146.3	210.1
1	7	271	508	175.8	391.5	567.5	23181.2	1502.2	11095.5	158.9	605.8
1	7	195	503	236.4	389.5	485.4	16827.3	6096.1	8300.3	189.5	5778.0
1	7	301	592	181.8	378.8	553.6	22568.7	8013.1	13222.2	889.6	1211.8
1	7	248	548	197.7	525.1	543.6	16512.2	1560.8	10651.3	175.3	508.9
1	7	242	532	208.7	386.0	508.1	16468.5	1355.0	8900.2	155.7	283.8



Bustard

## Sequence Files


s\_1\_0007\_seq.txt - WordPad

1	7	135	563	TTTGAACAAGCATATTGATAGCAGCAC
1	7	180	621	TGTTTTTTTTTTTTTTTTTGAGACAGAG
1	7	245	626	TTTGATCATGTTTTCTGCTGCTGAAGC
1	7	241	509	TCTCCTGCCTCAGCCTCCCAGTAGCT
1	7	214	595	TACAAAATCCCTGCCCATATGGAGCTT
1	7	155	544	TTATCTGCATCCGGTGCAAGTTTTAGC
1	7	301	507	TC CCTGCTTATTGACTCTTTTTTATTT
1	7	175	606	TTGGAATCGGGGTTAAAGGAAAGAGAT
1	7	242	522	TAACTAATATACAGGATATGTTCAAAA
1	7	196	522	TGTCACAGGAGGGAACAGCGCTGACAT
1	7	237	612	TTGCTGCAAGCTCAGAAGAACACTTTC
1	7	160	528	TCTGATTTTACACAGTAACAGAAAAC
1	7	164	543	TCTCAGAGAAACGTGCGTGATCCAGG
1	7	179	581	TTCTGAAATTACTGCCTGTCTATTGG
1	7	226	623	TATTACAGGCATGAGCCACTGCACCCA
1	7	139	583	TGTGGGTATGGGACACACAGGGAAGCT
1	7	220	618	TCGCAAGTTGTTTAAATAAGAGCAA
1	7	360	507	TTATTTGTGAGTAAATGTTTCCAATTA
1	7	334	512	TAGTTGGTTGCACCTAAATGGGAGATC
1	7	155	517	TCCACAAAAAGAAAAAAGAGAGAGA
1	7	343	541	TATGTTCCATGTGCTAATGAATAGAAT
1	7	241	608	TATTAGCCAGGTGTGGTGGTGTACACC
1	7	176	520	TTTTTTAGTAGAGATGGGATTTCAACA
1	7	371	592	TATTGCTATAGGAACAGCCAGTAGGGG
1	7	271	508	TCTCTGGGAAATATTAGCTTAGCCAGA
1	7	195	503	TACATGATGTGGGCCTGGTGATCTTG
1	7	242	532	TTCCTGCTGCTGCTGCTGCTGCTGCT

- Transforms intensities into base-calls

## From Reads to Aligned Sequences

## sequence files



Gerald

### Lane Results Summary

Line	Clusters	Area (a.u.)	% annually after 20 cycles	% PF Clusters	% Abiga (PF)	Area/alignment score (PF)	% Error Rate
1	36790 2495	1190 175	62.2 ± 0.8	49.3 ± 0.6	95.0 ± 0.2	1011.48 (8.7)	0.86 ± 0.2
2	36811 1862	1558 08	62.2 ± 0.8	47.2 ± 0.2	95.0 ± 0.4	2089.58 (8.62)	0.82 ± 0.0
3	37212 2322	1522 146	60.75 ± 0.6	49.75 ± 7.3	95.15 ± 0.4	3030.94 (8.68)	0.78 ± 0.0
4	36810 1768	1497.56	64.47 ± 2.0				
5	36801 1386	1427.92	60.5 ± 1.8				
6	37561 1362	1437.34	66.88 ± 0.87				
7	36802 1502	1524.30	66.16 ± 3.2				
8	36576 1719	1344 133	60.10 ± 0.7				

### Expanded Lane Summary



Lane Info		Plating Info	
Lane	Clonotype	% Plating	% Preplating
1	16799	0.650	0.630
2	18431	0.650	0.630
3	17212	0.650	0.630

**Lane 1**

Lane	Title	Clonotype	Av. Log10 CFU/50
1	00101	00000	1.00E+05
2	00102	00000	1.07E+05
3	00103	00000	1.22E+05
4	00000	00000	1.00E+05

**Layer 1**

Index	Title	Checksum	Age (in Days)
1	0011	00000	1000000
2	0012	00000	1000000
3	0013	00000	1000000
4	0014	00000	1000000



Export.txt






# CASAVA

## Consensus Assessment of Sequence and Variation

### Aligned Reads




Consensus assembly

BINS  
SNPs  
Counts



# Consensus Sequence





# Output Interpretation - html

● Project directory/html/Home.html



Human Genome


Statistics

Human Genome


## Human Genome

Run name	cycles columns	PF yield (whole FC)	Used Yield (Gb)	use lanes	Summ
<a href="#">080328_EAS114_0051_FC3011HAAXX_R1</a>	N/A	N/A	N/A	2, 3, 4, 6, 7, 8	✓
<a href="#">080328_EAS139_0029_FC3010UAAXX_R1</a>	N/A	N/A	N/A	1, 2, 3, 4, 6, 7, 8	✓
<a href="#">080328_EAS188_0025_FC20DRTAAXX_R1</a>	N/A	N/A	N/A	1, 2, 3, 4, 6, 7, 8	✓
<a href="#">080401_EAS192_0026_FC20DU9AAXX_R1</a>	N/A	N/A	N/A	1, 2, 3, 4, 6, 7, 8	✓
<a href="#">080401_EAS54_0088_FC20DT1AAXX_R1</a>	N/A	N/A	N/A	1, 2, 3, 4, 6, 7, 8	✓
<a href="#">080403_EAS114_0053_FC301CCAAXX_R1</a>	N/A	N/A	N/A	1, 2, 3, 4, 6, 7, 8	✓
<a href="#">080404_HWI-EAS6_0002_FC30178</a>	N/A	N/A	N/A	3, 4, 6, 7, 8	✓
<a href="#">080404_USI-EAS39_0005_FC301E9</a>	N/A	N/A	N/A	2, 6, 7, 8	✓
<a href="#">080407_SLXA-B7_0625_FC20DRKAAXX_R1</a>	N/A	N/A	N/A	1, 2, 3, 4, 5, 6, 7, 8	✓
<a href="#">080408_HWI-EAS300_0002_FC301FW</a>	N/A	N/A	N/A	1, 2, 6, 7, 8	✓
<a href="#">080409_SLXA-B7_0619_FC20DREAAXX_R1</a>	N/A	N/A	N/A	5, 6, 7, 8	✓
<a href="#">080409_SLXA-B7_0619_FC20DREAAXX_R1</a>	N/A	N/A	N/A	1, 2, 3, 4	✓
<a href="#">080411_USI-EAS39_0002_FC301KN_PE</a>	N/A	N/A	N/A	1, 2, 3, 4, 6, 7, 8	✓
<a href="#">080414_EAS192_0028_FC20DRMAAXX_R1</a>	N/A	N/A	N/A	1, 2, 3, 4, 5, 6, 7, 8	✓
<a href="#">080414_HWI-EAS218_0002_FC310MH_R1-R2</a>	N/A	N/A	N/A	1, 2, 3, 6, 7	✓
<a href="#">080414_HWI-EAS220_0003_FC310FK_R1-R2</a>	N/A	N/A	N/A	1, 2, 3, 6, 7, 8	✓
<a href="#">080415_HWI-EAS26_0003_FC301LW_R1-R2</a>	N/A	N/A	N/A	1, 3, 6, 7, 8	✓
<a href="#">080418_HWI-EAS6_0002_FC301KY_PE</a>	N/A	N/A	N/A	4, 6, 7, 8	✓
<a href="#">080418_USI-EAS39_0001_FC301KM_PE</a>	N/A	N/A	N/A	1, 2, 3, 4, 6, 7, 8	✓
<a href="#">080423_HWI-EAS218_0001_FC30307</a>	N/A	N/A	N/A	1, 3, 4, 6, 7, 8	✓
<a href="#">080423_HWI-EAS220_0001_FC30310</a>	N/A	N/A	N/A	1, 3, 4, 6, 7, 8	✓

Depth for all chromosomes



Number of SNP's for all chromosomes




# Overview of CASAVA Outputs for GenomeStudio™

- Output utilized for DNA experiments
  - Sorted.txt files binned by chromosome and subdivided into 10 megabase bins
  - SNP.txt (1 per chromosome)
  - Run\_summary.xml
  - Run.conf file
  - Project.conf file
- Output utilized for RNA experiments
  - Sorted.txt files binned by chromosome and subdivided into 10 megabase bins
  - SNP.txt (1 per chromosome)
  - Exon counts file (1 per chromosome)
  - Gene counts file (1 per chromosome)
  - Splice Junction counts file (1 per chromosome)
  - Run\_summary.xml
  - Run.conf file
  - Project.conf file


# DNA Sequencing Module

- Direct import of data from Pipeline/CASAVA
- Visualization of SNPs
- Browsing of coverage and consensus reads
- Export of SNP tables




## Exon Counts from RNA Data

## Browser view of exon counts



# A New Look at Alternative Splicing



# Leveraging the GA Informatics Community

## *De Novo Assembly*

- Velvet – De novo assembly of short reads
  - Daniel Zerbino and Ewan Birney, EMBL-EBI
  - <http://www.ebi.ac.uk/~zerbino/velvet/>
- SSAKE – Assembly of short reads
  - Group: Rene Warren, et al; British Columbia
  - <http://bioinformatics.oxfordjournals.org/cgi/content/full/23/4/500>
- Euler SR – Genomic Assembly
  - Group: Pavel Pevzner, Mark Chaisson; UC San Diego
  - <http://nbcrc.sdsc.edu/euler/>

# Rapidly Expanding Choice of Open Source Tools

## *Genomic Alignment Browsers*

- Gbrowse – Genomic Browsing
  - Generic Model Organism Database Project
  - <http://www.gmod.org/wlk/index.pho/Gbrowse>
- UCSC Browser – Genome browsing and comprehensive annotation
  - Generic Model Organism Database Project
  - <http://www.genome.ucsc.edu/goldenPath/help/customTrack.html>
- Anno-J – Genome Annotation and Visualization
  - Computational Systems Biology Center of Excellence
  - [http://www.anoj.org/csb\\_index.shtml](http://www.anoj.org/csb_index.shtml)



# Leveraging the GA Informatics Community

## *Alignment and Polymorphism Detection*

- MAQ – Mapping and Assembly with Quality
  - Heng Li, Sanger Centre
  - <http://maq.sourceforge.net/maq-man.shtml>
- SOAP – Short Oligonucleotide Alignment Program
  - Ruiqiang Li, Beijing Genomics Institute
  - <http://soap.genomics.org.cn/>
- Consed – Alignment and Polymorphism Detection
  - Green Lab, U. Washington (commercial offering)
  - <http://bozeman.mbt.washington.edu/consed/consed.html>

# Rapidly Expanding Choice of Open Source Tools

## ChIP Sequencing

- ChIP-Seq Peak Finder
  - Barbara Wold, Cal Tech and Rick Meyers, Stanford University
  - <http://woldlab.caltech.edu/html/software/>

## Digital Gene Expression

- Comparative Count Display
  - Alex Lash, NIH
  - <ftp://ftp.ncbi.nlm.nih.gov/pub/sage/obsolete/bin/ccd/>
- SAGE DGED Tool
  - Cancer Genome Anatomy Project
  - [http://cgap.nci.nih.gov/SAGE/SDGED\\_Wizard?METHOD=SS10,LS108ORG=Hs](http://cgap.nci.nih.gov/SAGE/SDGED_Wizard?METHOD=SS10,LS108ORG=Hs)

# Genome Studio: Data Analysis Platform for Many Applications

A single workbench with a growing number of modules

- Sequencing
- Genotyping (GT)
- Gene Expression (GX)
- Regulation (M)
- ChIP Sequencing (CS)
- We are changing the name to reflect a new platform →
  - microarray + sequencing



# Connecting with the larger informatics universe



- Illumina has a 3<sup>rd</sup> party partnership program – illumina•Connect – launched last year designed to “increase software/hardware ecosystem connecting Illumina to researchers in genomics, genetics and sequencing communities”
  - ~30 vendors and academic partners in the program
  - <http://www.illumina.com/pagesnrn.ilmn?ID=229>
  - Integrated tools connecting GenomeStudio to several 3<sup>rd</sup> party apps for microarray and sequencing data analysis

INGENUITY  
SYSTEMS

InforSense  
Embedding Intelligence Throughout the Enterprise

GENESIFTER  
.....  
Understand the Biology.

jmp  
Statistical Discovery™ From SAS.

infoQuant<sup>Q10</sup>  
Informatics with Intelligence

PROGENY

ROSETTA  
BIOSOFTWARE

Sapio  
Sciences

Partek  
turning data into discovery™

D N A S T A R

geospiza

BC  
Biocomputing  
Platforms Ltd

GLS  
GenoLogics

Genedata  
solutions in-silico

BIOTEAM  
Enabling Science

GenomeQuest

SNP  
quanti

gene

CLC bio


synaMATIX

illumina®

# Delivering on Roadmap Milestones

15x increase in 2008

4-5x increase in 2009





# Sequencing Enabled iScan

What is needed?

How does it work?

What can it do?

# iScan Sequencing Module technical Update



- Add on to enable sequencing on the iScan
  - Uses lasers and optics of iScan for imaging of flow cell
  - Fluidics module holding sequencing reagents, pumps and reagent chiller
- Sequencing specs
  - Throughput: ~ 0.5 GB per day, up to 5 GB per run
  - Data density: ~ 32 M clusters per run, 60M Reads
  - Flexible read length: from short single reads to 2x75 bb reads
  - 8 lanes, 1-12 samples per lane
- Applications supported
  - smallRNA, ChIP-Seq & mRNA seq
  - Targeted re-sequencing & re-sequencing of small genomes
  - High Density GWAS, Medium and HD Custom GTP (Array Applications)
  - WG Gene Expression Profiling, Human, Mouse and Rat (Array Applications)



# Thank You

Marco Cappelletti

[mcappelletti@illumina.com](mailto:mcappelletti@illumina.com)